


Encore Manual



enovate
by ergotron®

Contents

COMPANION GUIDES.....	3
WELCOME.....	4
PRODUCT MODELS COVERED BY THIS DOCUMENT	4
SYMBOLS.....	6
WARNINGS	7
WARNINGS	8
STATEMENT OF USE	9
SAFETY GUIDELINES.....	9
TRANSPORT / STORAGE ENVIRONMENT	10
STORAGE REQUIREMENTS	10
OPERATING ENVIRONMENT	10
WARRANTY	10
COMPLIANCE	11
EMI (ELECTROMAGNETIC INTERFERENCE).....	11
UNPACKING	13
OVERVIEW	15
WEIGHTS.....	17
CPU COMPARTMENT SIZE	17
I/O PORTS	17
SETUP.....	18
SIGHTLINE MONITOR ARM INSTALLATION.....	18
SIGHTLINE MONITOR ARM TENSION ADJUSTMENT	19
CPU INSTALLATION.....	20
VOLTAGE SETTINGS	20
KEYBOARD INSTALLATION.....	21
MOUSE INSTALLATION	21
OPERATIONS.....	22
SIGHTLINE MONITOR ARM MOVEMENT	24
HEIGHT ADJUSTMENT LEVER.....	25
RESTART PC.....	25
MOBIUSPOWER®	26
BATTERY PACK FUEL GAUGE	26
TIME REMAINING ON BATTERY PACK - LCD SCREEN SEQUENCE.....	27
MOBIUSPOWER® BATTERY PACK REMOVAL - LCD SCREEN SEQUENCE	27
MOBIUSPOWER® CHARGING STATION.....	28
SWAPPING THE MOBIUSPOWER® BATTERY	29
ON-BOARD CHARGING: (Encore MobiusPower® Plus Model only)	29
BATTERY DISPOSAL	30
POWER SPECIFICATIONS.....	31
ENCORE	31
CHARGING STATIONS	31
MOBIUSPOWER® BATTERY PACKS.....	31
MAIN POWER CONTROLLER ENCORE.....	32
MOBIUSPOWER® BATTERY HOLSTER.....	32
SYSTEM ELECTROMAGNETIC EMISSIONS	33
CLEANING	39
UI TOUCH SCREEN	39
WORK SURFACE.....	39
GENERAL WORKSTATION CLEANING INSTRUCTIONS.....	39
TESTED AND APPROVED CLEANING AGENTS	39
CLEANING PRECAUTIONS	40
MAINTENANCE	41
INSPECTING CORDS AND CONNECTORS.....	41
SAFETY CHECKS BEFORE ENERGIZING EQUIPMENT	41



COMPANION GUIDES

Encore Specification Sheet
Workstation Cleaning Instructions
Software User Guide
<https://software.ergotron.com>

WELCOME

Congratulations on your purchase of the Encore Mobile EHR Workstation. Please do not use your workstation until you have read this manual in its entirety as it contains important safety and use information. The Encore workstation should only be used as outlined in this manual. Be sure to keep this manual in a safe place for future reference. If at any time you have questions or concerns about the contents of this manual or the use of your Encore workstation, please contact Ergotron customer care.


This manual, product informational and safety labels related to the Encore workstation frequently employ the use of symbols with or without accompanying text. Health products regulatory agencies require the use of symbols, often in place of textual statements, to enhance the legibility of labeling and thus improve the conspicuousness of required information such as important safety information. Your attention to the presence and content of symbols used in this manual will help to ensure the safe use of your Encore workstation. The Symbols and Warnings chart depicts the symbols used and provides a definition for each symbol found in this manual and in the labels and labeling materials of the Encore workstation

PRODUCT MODELS COVERED BY THIS DOCUMENT





Model

- EC-MDLD
- EC-MDLN
- EC-MDSD
- EC-MDSN
- EC-MDAD
- EC-MDAN
- EC-MALD
- EC-MALN
- EC-MASD
- EC-MASN
- EC-MAAD
- EC-MAAN

Model information is located next to REF on the product serial label.



Encore





















WARNINGS

SYMBOLS

The term “symbols” refers to the use of graphical symbols without equivalent accompanying text. Symbols are used by medical device manufacturers to create uniform labels and labeling for the United States, European Union, and any other countries that may permit their use in medical products.

The use of symbols on the product labels are intended to conform to international consensus standards. The following chart depicts the symbols used and provides a definition for each symbol found on the labels.

	Alternating Current		
	Direct current		
	Caution In case of application as a safety sign, the rules according to ISO 3864-1 are to be adhered to. See safety sign ISO 7010- W001		
	Operating Instructions Refer to instruction manual/booklet		
	Dangerous Voltage		
	Caution, risk of electric shock		
	Date of Manufacture		Serial Number
	Manufacturer		Reference number of the device
	WEEE Symbol		Specific configuration of the device
	General Symbol for Recycle		Temperature range
	Arch Flash		Relative humidity
	No Pushing, to Prohibit Against an Object		Atmospheric Pressure

WARNINGS



Important Warnings



Electrical Shock Warning

The above symbols represent safety warnings that require significant attention when seen on the Encore Workstation or in the user manual. Failure to do so could result in minor injury, major injury, or even death.



NOT SUITABLE FOR USE IN AN OXYGEN RICH ENVIRONMENT! Do not use in the presence of an anesthetic mixture with air or with oxygen or nitrous oxide.



USE A NONFLAMMABLE CLEANER WHEN CLEANING THE UNIT! Failure to do so can result in death, explosion, and/or fire.



DO NOT LEAVE THE UNIT UNATTENDED AROUND CHILDREN! Failure to do so can result in injury, and/or death.



CAUTION: MAIN BATTERY IS REMOVABLE BY HANDLE. You must use proper lifting techniques. Failure to do so can result in injury.



AVOID USING AN EXTENSION CORD WITH THE UNIT! If an operational error occurs, the plug should be immediately removed from the socket.



THIS CHARGER IS DESIGNED FOR USE WITH LITHIUM BATTERIES!



For safety reasons, this charger must be used only for MobiusPower® batteries



DO NOT ATTEMPT TO SERVICE OR REPLACE ANY PART OF THE ENCORE WORKSTATION unless directed to do so through Ergotron approved documentation (i.e., this User Manual or other instructions). Only Ergotron or an Ergotron-certified entity may service or replace the cart components. If any component on the cart is missing or damaged, the cart must not be used. Contact Ergotron immediately to request service.



MAINTENANCE During preventive maintenance workstation should be turned off by removing main battery or disconnecting holster charger.



Tipping Hazard: Transport cart in lowered position using label placed in the lower column as a guide.

WARNINGS

DO NOT TRANSPORT THE WORKSTATION UP OR DOWN STAIRS! Workstation must be in lowered position and all drawers must be closed during transport.



DO NOT OPEN THE POWER SYSTEM! Unauthorized personnel opening the power system may cause injury and/or death. If the unit is not working properly, please contact Ergotron at 1-800-888-8458



DO NOT OPEN THE POWER SYSTEM! Unauthorized personnel opening the power system may cause injury and/or death. If the unit is not working properly, please contact Ergotron at 1-800-888-8458



DO NOT USE THE UNIT IN/NEAR WATER OR OTHER LIQUIDS! If the unit becomes wet, unplug it immediately, wipe away any excess liquid and allow it to dry before use. Failure to do so may cause electric shock, damage to the unit, voiding the warranty, injury or death.



ALWAYS KEEP THE UNIT WELL VENTILATED! Do not block ventilation airways or insert items into the ventilation slots. Failure to do so can cause the power system to overheat and possibly cause fire, explosion, and/or death.



THE CHARGER CONTAINS DANGEROUS VOLTAGES AND THE COVER SHOULD NOT BE REMOVED. All service or maintenance work should be carried out by qualified service personnel by contacting Ergotron at 1-800-888-8458



GROUNDING Connect the workstation or bay charger to an equivalent receptacle marked "Hospital Only" or "Hospital Grade" to ensure ground



WARNING: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth



WARNING: Do not push the workstation from the side. Always use the handle to move or adjust the workstation position.

STATEMENT OF USE

The Encore Mobile EHR Workstation facilitates mobile computing in healthcare environments to aid in clinical documentation. Designed to set a new standard in quality, ease of use, and customer satisfaction, our goal is to provide a product that is built to exacting standards and is ready for years of durable service in a healthcare environment. This product is designed to be safely used within general patient areas and is meant to aid in the entering or retrieving of clinical data. It complies with IEC 60601-1 electromagnetic leakage and safety requirements if used in accordance with the boundaries and suggestions of this manual.

SAFETY GUIDELINES

- Workstation is designed to ensure both the highest level of product quality and safety for the user. To maintain both quality and safety, follow the guidelines and instructions in this manual.
- Use the Workstation system only as intended.
- Do not place the Workstation near a window. Exposing the Workstation to rain, water, moisture, or constant, direct sunlight can severely damage it.
- The Workstation has no user-serviceable, internal parts. To maintain your warranty, refer all servicing to Ergotron qualified personnel.
- Do not cover or obstruct any venting holes on the Workstation or the MOBIUSPOWER® Charging Stations.
- Store the Workstation within 10 to 30 degrees Celsius (50 to 86 degrees Fahrenheit) for optimum backup battery life. Storing the Workstation outside the temperature range could result in premature backup battery failure.
- Use and maintain the cord set provided with the equipment.
- Position Workstation in a manner that does not obstruct or make it difficult to disconnect from external power source.
- If any cord or cable is frayed or damaged, replace it immediately with another of the same type and rating as supplied by Ergotron.
- To clean the exterior of the Workstation, follow the IEC 60601-1 standard for use in a hospital environment. See “Maintenance” for more information.
- Before cleaning a MOBIUSPOWER® Charging Station, disconnect the enclosure from its power source.
- CAUTION! DO NOT ship individual MOBIUSPOWER® Battery Packs by air. Certain restrictions apply. Contact Ergotron for shipping instruction.

CLASSIFICATIONS

- Charger: Class 1
- Main Controller: Class 1, Internally Powered
- Degree of Protection against Harmful Ingress of Water, IPX0
- EQUIPMENT not suitable for use in the presence of a FLAMMABLE ANAESTHETIC MIXTURE WITH AIR or WITH OXYGEN OR NITROUS OXIDE

TRANSPORT / STORAGE ENVIRONMENT

Care should be taken to transport and store this system within the following:

- Ambient Temperature Range: -20°C to +50°C (-4°F to 122°F) Storage
- Atmospheric Pressure: 50 kPa to 106 kPa
- Relative Humidity: 5% to 85% non-condensing

STORAGE REQUIREMENTS

The lithium ion Battery Pack must not be fully charged prior to storing. 50% state of charge recommended.

OPERATING ENVIRONMENT

- Ambient Temperature Range: 0°C to +40°C (32°F to 104°F)
- Atmospheric Pressure: 50 kPa to 106 kPa
- Relative Humidity: 5% – 85% non-condensing

WARRANTY

Product Warranty:

2 year limited warranty on Battery Packs

5 years on Workstation and electronic components

Contact Ergotron directly for full warranty details

POWERING DOWN

To power off the Encore workstation:

1. Unplug power cord from power source (if applicable)
2. Remove the swappable battery pack
3. Press and hold the power button for >10 seconds

COMPLIANCE

Workstation

- IEC 60601-1

Bay Chargers

- UL 1012
- CAN/CSA-E60335
- IEC 60335

Mobius Battery

- IEC 62133
- UN 38.3

EMI (ELECTROMAGNETIC INTERFERENCE)

Portable and mobile RF communications equipment can affect Medical Electrical Equipment. The use of accessories, transducers, and cables other than those specified by the manufacturer, may result in increased Emissions or decreased Immunity of the System. The System should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the System should be observed to verify normal operation in the configuration in which it will be used.


These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instruction manual, may cause interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio, television or Medical Electrical Equipment reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference using one or more of the following measures: reorient or relocate the receiving antenna; increase the separation between the equipment and the receiver; connect the equipment into an outlet on a circuit different from that which the receiver is connected; consult the dealer or an experienced radio/television technician for help. The user must use shielded cables and connectors with this product.

Any changes or modifications to this product not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation. MOBIUS POWER meets or exceeds FCC Class A limits for EMI

UNPACKING


Promptly unpack your products to check for completeness and damage caused by shipping. Immediately after receiving your products, ALL batteries must be fully charged to ensure the duration of their warranty.

NOTE: All Pictures shown are for illustration purpose only. Actual product may vary due to product enhancement.




STEP 1

Use scissors or a utility knife to cut and remove the two outer straps.

**STEP 2**

Remove cardboard lid, cardboard wall and any packing spacers.



NOTE: STEP 3 requires two people.


STEP 3

Lift cart (with foam caster braces still attached) out of the cardboard base and remove foam caster braces.



STEP 4


Open tray and check for integration kit and locate all optional accessories

OVERVIEW



*OPTIONAL FEATURE BASED ON MODEL



**WEIGHTS**

Total Mass = 70.3 kg

*Includes battery, power and maximum accessory limits.

Monitor Pole = 7.3 kg (16 lbs.)

Monitor Arm = 4.1 kg (9 lbs.)

CPU Compartment = 4.5 kg (10 lbs.)

Work Surface = 4.5 kg (10 lbs.)

Keyboard Tray = (3 lbs.)

Column Accessories = 4.5 kg (10 lbs.)

Column Drawer = 9.1 kg (20 lbs.)

CPU COMPARTMENT SIZE

Width: 406.4mm (16")

Depth: 330.2mm (13")


Height: 57.2mm (2.25")

I/O PORTS:


Channel 1 - PC

Channel 2 - Monitor

Channel 3 - Peripheral




SETUP



SIGHTLINE MONITOR ARM INSTALLATION

STEP 1


Align screws with the four holes in the monitor pole VESA plate. Hold the monitor in place and tighten the screws.




SIGHTLINE MONITOR ARM TENSION ADJUSTMENT

STEP 2

Use an Allen Wrench to adjust the SightLine Monitor Arm tension, turn clockwise to loosen, turn counter clockwise to tighten.

**STEP 3**

Connect the monitor's power cable and the video cable. Route each cable from the monitor through the cable management cover located on the front of the monitor arm. Leave enough slack above the cable management cover to allow the monitor to extend fully upward.

**STEP 4**

Continue routing each cable behind the arm and through the rear cable port at the base of the monitor arm and plug into computer port.


CPU INSTALLATION


Tilt Monitor back, unlock worksurface, and open worksurface. Set the CPU onto the platform beneath the work surface. NOTE: Take care in CPU placement; consider orientation based on direction of cable inputs and sizes INCLUDING monitor input, power input, mouse, keyboard and any peripherals. Make cable connections within the tray.

VOLTAGE SETTINGS can be changed using EM Control or the Rhythm Unplugged app. For additional information please visit: <https://software.ergotron.com>

I/O PORTS

- Channel 1 - PC
- Channel 2 - Monitor
- Channel 3 - Peripheral





KEYBOARD AND MOUSE INSTALLATION

STEP 1

Set the keyboard on the keyboard tray

STEP 2


Pull the keyboard and mouse tray to full extension. This will simulate the full length requirement of the mouse cable length.

STEP 3

Route cables under and through the tray routing ports, plug into appropriate location on CPU/laptop and zip tie into place.

STEP 4

Replace the work surface and cable cover at the base of the monitor arm, close and lock.



STEP 5

Stow mouse in mouse holder

OPERATIONS

POWERING ON


To power on the Encore, insert a fully charged MobiusPower® Battery Pack or plug into an outlet.

As the Encore is powering on, you will see a series of screen displays. For full details please visit:
<https://software.ergotron.com>




SWAPPING THE MOBIUSPOWER® BATTERY

Remove MobiusPower® battery pack from holster, swap with a fully charged battery from the charging station.




SIGHTLINE MONITOR ARM MOVEMENT

Upright Position




Extended Position




Reclined Position

For safe travel, lower the Sightline Monitor Arm to the reclined position.



Monitor rotates 80 degrees.



HEIGHT ADJUSTMENT LEVER

To raise or lower the work surface:

STEP 1


Position your hands to hold the work surface steady and grip the height adjustment lever located on the right side of the Encore,

STEP 2

Pull the height adjustment lever upward to raise or lower to desired position.

STEP 3

Release the height adjustment lever once the desired position is reached.









Up/Down Height
Adjustment Lever

MOBIUSPOWER®

BATTERY PACK FUEL GAUGE

- Each Battery Pack has a “fuel gauge” to indicate charge level/capacity.
- You may check a Battery Pack’s charge capacity at any time.
- To check the remaining charge on a Battery Pack, press and release the white circle on the Battery Pack’s fuel gauge display.
- The charge level is indicated in 25% increments by the number of status lights appearing between the E-Empty and F-Full.
- After 3 seconds, the status lights will turn OFF.




Charge Capacity		LED 1	LED 2	LED 3	LED 4
	0% - 10%	FLASH	OFF	OFF	OFF
	11% - 25%	ON	OFF	OFF	OFF
	26% - 50%	ON	ON	OFF	OFF
	51% - 75%	ON	ON	ON	OFF
	76% - 100%	ON	ON	ON	ON

TIME REMAINING ON BATTERY PACK - UI SCREEN

The workstation is calculating the run time based on the current being drawn by the system.



Once the time is calculated, the UI will show a continuous display of 'Time Remaining' in hours and minutes. The time will vary depending on usage amount.




When the battery pack is low, UI will show notification to swap or charge battery soon.


MOBIUSPOWER® CHARGING STATION

Includes 2-Battery or 4-Battery Pack Charging Station. Once the Charging Station is placed and plugged into an outlet, slide a Battery Pack into any open slot making sure the FLASHING GREEN indicator light turns on. The charge cycle is complete when the indicator light is illuminated SOLID GREEN.

MobiusPower® Dual Pack Charging Station



Charge Indicator Lights



FULLY CHARGED

Solid green light



CHARGING

Blinking green light




ERROR

Solid or blinking yellow light

MOBIUSPOWER® HOLSTER CHARGER

To charge the MobiusPower® battery pack in the holster, plug the power cord into an AC electrical outlet. The UI message center will read “AC Plugged In Charging”.

The charge status can be checked by unplugging the AC power cord from the wall outlet and reading the time remaining in hours and minutes on the UI message center or depress the charge level indicator on the MobiusPower® battery pack.




MOBIUSPOWER® BATTERY PACK SWAPPING OR REMOVAL - UI SCREEN

When the battery pack is removed from the workstation holster, the workstation will begin a 2 minute countdown timer. The workstation will temporarily draw power from the DC outputs and internal backup battery.

The workstation will emit an audible alert every 5 seconds during the 2 minutes.

If the Battery Pack charge capacity is depleted or no fully charged battery is attached, the UI will turn off. At this point, MobiusPower®, your computer, and any other devices requiring power will shut down until the battery pack is replaced with a Fully Charged Battery, OR the cart is plugged into an AC outlet. You may then turn on your computer and other electronic devices.

ON-BOARD CHARGING: (Encore MobiusPower® Plus Model only)

The Encore MobiusPower® Plus can charge the MobiusPower® battery pack while the battery pack is in the holster. To begin charging your Encore MobiusPower® Plus, plug the 3-prong power cable into a standard 110-120VAC grounded AC power outlet.

The UI will show charging when the AC cord is connected. Please unplug the workstation from the AC Power outlet prior to moving it.

BATTERY DISPOSAL

Industrial batteries contain materials which are considered 'hazardous substances'. If batteries are improperly disposed of, for example, thrown in the trash or illegally dumped, these substances can eventually leak out and contaminate the surrounding soil and groundwater supply.

It is a violation of federal law to improperly dispose of batteries once they can no longer be used. Once a battery is purchased, full liability and responsibility lies on the owner to dispose of the battery. The law says that responsibility is still on the owner if the battery is disposed of improperly by dumping in a landfill, or shipping to a scrap dealer who does not handle it properly and in which environmental damage occurs.

It is illegal to dispose of batteries in any way other than 'thermal recovery' or recycling of the hazardous substances in batteries according to the Environmental Protection Agency (EPA).

The Department of Transportation (DOT) has strict guidelines for the shipping of hazardous materials, which result in large fines if they are not followed.

Check with your local ordinances for proper disposal procedures. SDS sheets are available on-line at <http://www.ergotron.com>. Batteries are consumable goods.

Recycling MOBIUSPOWER® Batteries - Rechargeable batteries power everything from portable phones and cordless phones to laptops and PDAs. Batteries are a unique product comprised of heavy metals which include nickel cadmium, alkaline, mercury, nickel metal hydride and lead acid, which can threaten our environment if not properly discarded or handled. There are many ways to properly dispose of batteries, most of which depends on the type of battery you have. Lithium-ion rechargeable batteries are less toxic than many others, but it is still recommended that they be recycled.

The Rechargeable Battery Recycling Corporation (RBRC) manages a collection and recycling program for rechargeable batteries. The RBRC accepts rechargeable nicad, NiMH, Lithium-Ion, and small (under 2-lb) sealed lead-acid batteries. The RBRC has made it very convenient for you to recycle your rechargeables by getting retail stores like Home Depot, Target, Wal-Mart, and others to serve as collection points. And it's all free of charge. To find a drop-off point near you in the US or Canada use the Call2Recycle Locator tool at <http://www.call2recycle.org/> (U.S.); <http://www.call2recycle.ca/> (Canada).

Remember: Batteries enable our mobility, so it's likely our society will be using a lot more batteries in the future. But to ensure that we're not slowly poisoning our highly mobile selves, it's important that we do a good job of recycling batteries.



Equipment & Accessories Disposal

1. Please dispose of all batteries in accordance with local law.
2. All Electronics should be recycled through an electronics recycler.
3. Remaining plastics and metals can be recycled through a commercial recycler.

POWER SPECIFICATIONS

ENCORE

Model	EC-MA	EC-MD
Input Voltage	120 VAC	10.8 VDC
Input Current	4.2 A max	26.1 Ah max
Input Frequency	60 Hz	DC
Output Voltage		
Channel 1	5 – 24 VDC	5 – 24 VDC
Channel 2	5 – 24 VDC	5 – 24 VDC
Channel 3	5 – 24 VDC	5 – 24 VDC
Output Current	5 A max, each channel	5 A max, each channel
USB Charging Port	5V, 2.4A	5V, 2.4A

CHARGING STATIONS

Charging Specs	Dual Charging Station	Quad Charging Station
Inlet	IEC 320 C14	IEC 320 C14
Protection	10 A, 250 VAC Fuse	15A BREAKER, 250VAC
Input Voltage	120V +/- 10%	120V +/- 10%
Input Frequency	60Hz	60Hz
Input Current	3.6 Amps max	7.2 Amps max
Output Voltage	12.6 VDC	12.6 VDC
Output Current	8.5 A per bay (max)	8.5 A per bay (max)
Output Power	115 Watts (max) per bay	115 Watts (max) per bay
Measurements	32 cm (12.6") H x 43.2 cm (17") W x 14.7 cm (5.8") D	32 cm (12.6") H x 68.6 cm (27") W x 14.7 cm (5.8") D

MOBIUSPOWER® BATTERY PACKS

Battery Specs	Two Year Battery
Chemistry	Lithium Ion
Capacity	26.1 AH
Voltage	10.8 (nominal)
Charge Time	< 4 hrs.
Measurements	31 cm (12.2") H x 12.7 cm (5") W x 8 cm (3.14") D
Weight	2.5 kg (5.5 lbs.)

MAIN POWER CONTROLLER ENCORE

Input Voltage	10.8 VDC (nominal)
Input Voltage Range	9.0 VDC – 12.6 VDC
Input Current	14.5 A (max)
Output Voltage	3 Channels variable 5 VDC – 24 VDC
Output Power	3 Channels 120 W (max) each
Output Power Total	161 W (nominal)
Output Protection	Yes
Overload Recovery	Yes
Overload Protection	Yes

MOBIUSPOWER® BATTERY HOLSTER

Holster Assembly	Encore
Input Voltage	8.5 VDC ~ 15 VDC
Input Frequency	DC
Input Current	30 A (max)

SYSTEM ELECTROMAGNETIC EMISSIONS

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.37

WARNING: Use of unapproved ACCESSORIES may result in degradation.

Guidance and Manufacturer's Declaration – Cables, Transducers and Accessories			
The listed cables, transducers and accessories have been determined by Ergotron to be compliant with the emissions and immunity requirements of IEC 60601-1-2: 2014.			
Part No.	Description	Max Possible Length	Shielded (Y/N)
P0000937	CABLE, IEC AC INLET	3048 mm (10')	N
ENDC0001L	ECOFLEX R7 DC PWR CABLE 50 IN	1270 mm (50")	N
ENDC0002L	ECOFLEX DC PWR CBL 2.5X5.5X9 50 IN	1270 mm (50")	N
ENDC0003L	ECOFLEX DC PWR CBL 2.5X5.5X14 50IN	1270 mm (50")	N
ENDC0004L	ECOFLEX DC PWR CBL 2.1 X 5.5 50IN	1270 mm (50")	N
ENDC0005L	ECOFLEX DC PWR CABLE MP205 50IN	1270 mm (50")	N
ENDC0006L	ECOFLEX DC PWR CBL 2.1X5.5X14 50IN	1270 mm (50")	N
ENDC0007L	ECOFLEX DC PWR CBL HP PC 150W or less 50IN	1270 mm (50")	N
ENDC0009L	CABLE, DC POWER 2.1X5.5 X 14 RP 50IN EF	1270 mm (50")	N
ENDC0011L	ASM, DC PWR HP PC 150W 4.5mm 50", ECOFLX	1270 mm (50")	N
ENDC0012L	CABLE, 2C ULTRFT 50IN 3X6.5X9.5 OVRMLD	1270 mm (50")	N
ENDC0013L	CABLE, LENOVO 65W USB-C ADAPTER, ECOFLEX	1270 mm (50")	N
ENDC0014	ECOFLEX DC POWER SPLITTER	330.2 mm (13")	N
ENDC0016	CABLE, DC PWR, LG 24IN IGEL AIO, ECOFLEX	1270 mm (50")	N
ENDC0017	CABLE, WELCH ALLYN CONNEX SPOT DC CABLE	1270 mm (50")	N
ENDC0018	CABLE, 90W USB-C ADAPTER, ECOFLEX PWR	1270 mm (50")	N

TABLE 201 – GUIDANCE AND MANUFACTURER’S DECLARATION – ELECTROMAGNETIC EMISSIONS – FOR ALL EQUIPMENT AND SYSTEMS

Guidance and Manufacturer’s Declaration – Electromagnetic Emissions		
Encore is intended for use in the electromagnetic environment specified below. The customer or the user of Encore should assure that it is used in such an environment.		
Emission Test	compliance	Electromagnetic environment guidance.
RF emissions CISPR 11 EN 5505: 2009 + A1: 2010	Group 1	Encore uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonic emissions IEC 61000-3-2 EN610003-2: 2014	Class A	Encore is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Voltage fluctuations / flicker emissions IEC 61000-3-3 EN 61000-3-3: 2013	All Parameters	

TABLE 202 – GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY – FOR ALL EQUIPMENT AND SYSTEMS

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
Ecore is intended for use in the electromagnetic environment specified below. The customer or the user of Ecore should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2 EN 61000-4-2: 2009		±2,4, 6 & 8 kV Contact ±2, 4, 8 & 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/ burst IEC 61000-4-4 EN 61000-4-4: 2012		±2 kV for power supply lines ±1 kV for input/ output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5 EN 61000-4-5: 2006		±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 EN 61000-4-11: 2004	<5% Ut (>95% dip in Ut) for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec <5% Ut (>95% dip in Ut) for 1 cycle 40% Ut (>60% dip in Ut) for 5 cycle	<5% Ut (>95% dip in Ut) for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec <5% Ut (>95% dip in Ut) for 1 cycle 40% Ut (>60% dip in Ut) for 5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ecore requires continued operation during power mains interruptions, it is recommended that the Encore be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 EN 61000-4-8: 2010	30/ A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: Ut is the a.c. mains voltage prior to application of the test level.			

TABLE 204- GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY – FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE-SUPPORTING

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
Encore is intended for use in the electromagnetic environment specified below. The customer or the user of Encore should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic Immunity Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Encore including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF	3 Vrms outside industrial, scientific and medical (ISM) and amateur radio bands. 6 Vrms in ISM and amateur radio bands	6 Vrms 150 kHz to 80 MHz	$d = [1.17]\sqrt{P}$
IEC 61000-4-6 EN 61000-1-6: 2014	150 kHz to 80 MHz		
Radiated RF	10 V/m	10 V/m	$d = [1.17]\sqrt{P}$...80 MHz to 800 MHz
IEC 61000-4-3 EN 61000-4-3: 2006+ A1: 2008 + A2:2010	80 MHz to 2.7 GHz	80 MHz to 2.7 GHz	$d = [2.33]\sqrt{P}$...800 MHz to 2.5 GHz
	27 V/m, 18 Hz PM 385 MHz 28 V/m, 50% 18 Hz PM 450 MHz	27 V/m, 18 Hz PM 385 MHz 28 V/m, 50% 18 Hz PM 450 MHz	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)

TABLE 204 CONTINUED- GUIDANCE AND MANUFACTURER'S DECLARATION –
ELECTROMAGNETIC IMMUNITY – FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE-
SUPPORTING


Guidance and Manufacturer’s Declaration – Electromagnetic Immunity			
Encore is intended for use in the electromagnetic environment specified below. The customer or the user of Encore should assure that it is used in such an environment			
Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Immunity Guidance
	9 V/m, 217 Hz PM 710 MHz	9 V/m, 217 Hz PM 710 MHz	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range. b) Interference may occur in the vicinity of equipment marked with the following symbol: 
	9 V/m, 217 Hz PM 745 MHz	9 V/m, 217 Hz PM 745 MHz	
	9 V/m, 217 Hz PM 780 MHz	9 V/m, 217 Hz PM 780 MHz	
	28 V/m, 18 Hz PM 780 MHz	28 V/m, 18 Hz PM 780 MHz	
	28 V/m, 18 Hz PM 810 MHz	28 V/m, 18 Hz PM 810 MHz	
	28 V/m, 18 Hz PM 870 MHz	28 V/m, 18 Hz PM 870 MHz	
	28 V/m, 18 Hz PM 930 MHz	28 V/m, 18 Hz PM 930 MHz	
	28 V/m, 217 Hz PM 1720 MHz	28 V/m, 217 Hz PM 1720 MHz	
	28 V/m, 217 Hz PM 1720 MHz	28 V/m, 217 Hz PM 1720 MHz	
	28 V/m, 217 Hz PM 1845 MHz	28 V/m, 217 Hz PM 1845 MHz	
	28 V/m, 217 Hz PM 1970 MHz	28 V/m, 217 Hz PM 1970 MHz	
	27 V/m, 217 Hz PM 2450 MHz	27 V/m, 217 Hz PM 2450 MHz	
	9 V/m, 217 Hz PM 5240 MHz	9 V/m, 217 Hz PM 5240 MHz	
	9 V/m, 217 Hz PM 5500 MHz	9 V/m, 217 Hz PM 5500 MHz	
9 V/m, 217 Hz PM 5785 MHz	9 V/m, 217 Hz PM 5785 MHz		
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects			
a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cord- less) telephones and land mobile radios, amateur radio, AM and FM radio broad- cast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the EcoFlex Power is used exceeds the applicable RF compliance level above, the EcoFlex Power should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating EcoFlex Power.			
b. Over the frequency range 150 kHz to 80 MHz, field strength should be less than [V1] V/m.			

TABLE 206 – RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE EQUIPMENT OR SYSTEM – FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the Encore			
Encore is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Encore can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Encore as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d = $[1.17] \sqrt{P}$	80 MHz to 800 MHz d = $[1.17] \sqrt{P}$	800 MHz to 2.5 GHz d = $[2.33] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.75
1	1.17	1.17	2.33
10	3.70	3.70	7.36
100	11.70	11.70	23.30
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

CLEANING

UI TOUCH SCREEN

The Encore UI Touch Screen may be cleaned with a soft cloth or tissue dampened with water (or a mild detergent-water solution). If a detergent solution is used, rinse with a clean tissue dampened with water only. **CAUTION:** The Encore UI is not watertight. Do not use abrasive wipes or tissues on the UI Touch Screen – abrasive wipes may scratch the window. Never use solvents on the window – solvents may damage the finish or Touch Screen.

WORK SURFACE

All configurations of Encore workstations utilize an external plastic covering that is designed to resist the effects of harsh chemicals in a health care environment.

GENERAL WORKSTATION CLEANING INSTRUCTIONS

Nonabrasive cleaners or mild cleaning solutions should be used. As a precaution, we recommend testing the suitability of a cleaning product by applying to an inconspicuous area, minimizing the time of exposure and the amount of cleaning agent used (dilute as recommended by the supplier) to prevent any damage to the surface.

TESTED AND APPROVED CLEANING AGENTS

- Quaternary ammonium compound
 - Bleach Solution (generic)
 - *Alcohol Solution (91% isopropyl: Avoid Use on Base)
 - Hydrogen Peroxide
 - Sani-Cloth® Plus Wipes
 - Super Sani-Cloth® Wipes
-
- Contamination by intensely colored substances such as: coffee, iodine, or dyes, must be removed immediately.
 - Power System should be inspected bi-annually to ensure vent holes and pan guard are free of dust and debris.
 - All paints and plastic cart components will withstand cleaning by most commonly used diluted non-abrasive solutions such as: quaternary ammonia compounds, or bleach.

CLEANING PRECAUTIONS

- Do not use steel wool or other abrasive cleaners, solvents, polishes, waxes or steam cleaning tools that will damage the surface finish.
- Do not use strong solvents such as trichloroethylene and acetone; these materials will damage the surface finish.
- Caution should be used in the application of cleaning agents containing Ethylene Glycol Monobutyl Ether including: CaviWipes™, 409® Cleaner, Windex® and ammonia cleaners which have been shown to cause increased brittleness in PC/ABS components over time.
- The use of Zymit® Enzyme Cleaner, Virex® 256, or Micro-90® is not recommended.
- To avoid risk of electric shock, do not expose electrical components to water, cleaning solutions or other potentially corrosive liquids or substances.
- Do not immerse the cart or cart components in liquid or allow liquids to flow into the cart. Wipe cleaner off surface immediately using a damp cloth. Thoroughly dry surface after cleaning.
- Do not use flammable cleaners on cart surfaces due to proximity of electrical power and equipment.

IMPORTANT INFORMATION

- The following recommendations should not be construed to guarantee infection control. The hospital infection control administrator should be consulted regarding cleaning procedures and processes.
- *In some cases, alcohol will strip paint stampings like custom branding and graphics

MAINTENANCE

No user serviceable parts, repairs and/or upgrades are to be performed by an authorized service provider only. Your device provides reliable and efficient operation with a minimum of care. Although specific maintenance is not required, the following periodic checks ensure dependable operation:

INSPECTING CORDS AND CONNECTORS

Equipment can be damaged during use and when it's moved from one location to another. Damage can render equipment unsafe to use. Device interface cables, both power and data, ship installed in the Encore workstation.

PERFORM THE FOLLOWING **SAFETY CHECKS BEFORE ENERGIZING EQUIPMENT**

- When properly installed the plug should fit snug and not wobble in the socket. Look for visible damage to equipment, wall disconnects (holster charger), enclosures and electrical devices
- **IMPORTANT!** Ensure power supply cord has proper grounding with a 3-prong plug
- Inspect cables for wear, look for defects and missing, bent or broken plugs
- Inspect external sheaths of power supply cords to ensure they are not cut, frayed, twisted, or damaged
- Ensure the inner cores of power supply cords are not exposed where they connect to equipment, plugs and/or sockets
- Ensure inner cores of power supply cords are not exposed or twisted along the cord's length
- Ensure there are no exposed, unprotected electrical conductors
- Ensure electrical tape or anything that can cover damage has not been applied to power supply cords and devices
- **NOTE:** If any of the above items exist, immediately tag the item and notify a supervisor or facilities management
- In Addition:
 - Switch off disconnects and equipment before connecting or disconnecting power supply cords
 - Keep power supply cords clear of walkways (holster charger)
 - Keep power supply cords away from heat, sharp edges and moving parts that can damage cord sheath and insulation
 - To prevent overheating, do not use extension cords with MobiusPower® Plus
 - Notify maintenance if a power supply cord feels more than comfortably warm
 - Never use adapters to convert 3-prong plugs into 2-prong plugs
 - Pull the plug, not the cord when unplugging equipment. Pulling the cord causes unnecessary wear and adversely affects wiring to the plug
 - Walk or look all around equipment prior to moving to ensure it has been completely disconnected from the wall outlet



www.ergotron.com | USA: 1-800-888-8458 | Europe: +31 (0)33-45 45 600 | China: 400-120-3051 | Japan: japansupport@ergotron.com

© 2025 Ergotron, Inc. All rights reserved. Enovate is a registered trademark of Ergotron Inc.

While Ergotron, Inc. makes every effort to provide accurate and complete information on the installation and use of its products, it will not be held liable for any editorial errors or omissions (including those made in the process of translation from English to another language), or for incidental, special or consequential damages of any nature resulting from furnishing this instruction and performance of equipment in connection with this instruction. Ergotron, Inc. reserves the right to make changes in the product design and/or product documentation without notification to its users. For the most current product information, or to know if this document is available in languages other than those herein, please contact Ergotron. No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise without the prior written consent of Ergotron, Inc., 1181 Trapp Road, Eagan, Minnesota, 55121, USA. Patents Pending and Patented U.S. & Foreign. Ergotron is a registered trademark of Ergotron, Inc.