



ENVOY 2.0 SPECIFICATIONS - MOBIUSPOWER®


DIMENSIONS



DIMENSIONS WITH MED STORAGE MODULE




INTERNAL DESK SPACE




Designed to accommodate a 15" - wide Laptop

AUTOTRAX™ CASTERS



ENVOY 2.0 SPECIFICATIONS - MOBIUSPOWER®

MONITOR




Monitor Rotation 80°

STRAIGHT ARM MONITOR MOVEMENT




45° Upward Tilt
15° Downward Tilt

SIGHTLINE ARM MONITOR MOVEMENT




Upright
Extended Forward
Reclined

KEYBOARD, WORKSURFACE, MOUSE SHELF



Maximum Functional Tilt 45°




Mouse House






Mouse Shelf Extends
Left or Right

LCD TOUCH SCREEN



MOBIUSPOWER



Envoy Mobile Workstation Weight: 95 lbs (powered)
Mobius Battery Weight: 5.5 lbs
Monitor Pole Weight Limit: 16lbs

www.enovatemedical.com